

January 16, 2026

Dr. Hyosub Kim, Assistant Professor
Department of Education | School of Kinesiology
The University of British Columbia
5955 University Boulevard,
Vancouver BC, V6T 1Z1

Via email: hyosub.kim@ubc.ca
Dear Dr. Kim/Hyosub,

It is an honour to welcome you back to the Kinarm User Community! We hope that your new Kinarm End-Point Lab will be an invaluable addition to your research program and provide you with many years of successful studies and publications. Now that we have completed the installation and training, there are a few things we would like to remind you of:

1. We are here to support you.

If you haven't already done so, please take a moment to "[register](#)" on the Learning & Support page of the Kinarm Website. With this login you will be able to access lots of information such as User Guides, and new releases of Dexterit-E and Dexterit-E Explorer. Anyone from your lab can request an account. We hope you and your students will take advantage of it to the fullest. We note you have created a single 'lab account' under the above email. This account will receive notifications of releases and product notices/advisories.

2. Control who uses your Kinarm Lab.

Remember that your Kinarm is a sophisticated (and expensive!) piece of instrumentation. All users in your lab should receive mandatory training before they begin to use the Kinarm. The Kinarm Operator Guide is a good place for a new user to start. While incidents of graduate students "breaking" a Lab are thankfully rare, it has happened. I would encourage you to exercise control over who has access to your lab. This control should also extend to persons outside of your lab, including the University IT department. We have seen a rise in institutions implementing cybersecurity measures that have adversely affected the Dexterit-E computer. We have seen a rise in institutions implementing cybersecurity measures that have adversely affected the Dexterit-E computer. **Please be vigilant in the early days of the system and ask support@kinarm.com before the IT department touches the Kinarm Lab.**

3. Keep the computers "clean".

The computers that form part of the Kinarm control system have been set up to ensure robust, reliable operation of your Lab. Changing settings on these computers can adversely affect the performance of the Kinarm Lab and may void the warranty.

BT088 Welcome Letter.Kim.16-Jan-26

Form: BT088 Welcome Letter.16-Jan-26

Please be sure that all users and administrators of your Kinarm Lab are aware that there are notable differences between a Dexter-E computer and a standard installation of Microsoft Windows (per Chapter 5, [Dexter-E User Guide](#)) such as:

- **Anti-Virus** - The Dexter-E computer comes with anti-virus software installed. The installed anti-virus software is known to work properly with Dexter-E and your Kinarm Lab. Many anti-virus software applications cause problems with Dexter-E. In particular, many anti-virus software applications interfere with the subject display's timing and, therefore, should not be used.
- **Backup** - See [Section: 5.3.6 Database Backup](#) for details on how to back up your data. If multiple Windows accounts are created, you may want to ensure that the read permission for the backup destination folders are set correctly to maintain data isolation and privacy.
- **Networking** - The Dexter-E computer should not be added to a Microsoft domain. Microsoft domains make significant changes to the permissions on the Dexter-E computer and can result in Dexter-E being unable to run.
- **Accounts** - You may create multiple Microsoft Windows accounts for users to provide data isolation between users and groups. Accounts do not require administrative privileges to run Dexter-E. Each new user account that you create will have a separate, default location where all relevant exam data are saved. See **Backup** above for issues related to multiple accounts and backups.

4. **Look after your new Kinarm Lab.** In the [Kinarm Lab Operator Guide](#) there is an Annual Maintenance Checklist. Please complete the checklist once a year to ensure your Lab continues to run "as good as new".

We look forward to supporting you and your research program in the years to come!

Kind Regards,



Digitally signed by Anne Vivian-Scott
Reason: I am the author of this document
Location: Kingston, Ontario, CANADA
Date: 2026.01.16 14:37:30 -05'00'

Anne Vivian-Scott, President & CEO
avivianscott@kinarm.com

BT088 Welcome Letter.Kim.16-Jan-26

Form: BT088 Welcome Letter.16-Jan-26

Revision History

Revision Date	Changes Made	Reason for Change	Approved By:
24-Aug-21	Made the standard letter a BT form. Incorporated changes from IT-136 and harmonized by Ch6 of Dex UG.	With growth and need to control interactions with customers per 72BT02	Anne Vivian-Scott
21-Sep-23	Updated links	new website & 3.10 release	Anne Vivian-Scott
16-Jan-26	Updated links; language re cyber protection	3.11 release & CAPA-970	Anne Vivian-Scott

BT088 Welcome Letter.Kim.16-Jan-26

Form: BT088 Welcome Letter.16-Jan-26